

Kentucky Method 64-254-083
Revised ~~1/6/03~~03/07/08
Supersedes 64-254-030
Dated ~~1/27/00~~01/06/03

CHEMICAL ANALYSIS OF CEMENT

1. SCOPE: This method is a modification of ASTM C-114 for the chemical testing of cement by x-ray spectroscopy.
2. APPARATUS –MATERIALS:
 - 2.1. Philips MagiX PRO Wavelength Dispersive x-ray fluorescence spectrometer.
 - 2.2. SuperQ software
 - 2.3. Philips Perl’X 3 fused bead machine
 - 2.4. Platinum dish and crucible set
 - ~~2.5. Spex 8000 mixer/mill~~
 - 2.56. Lithium Bromide (LiBr): 10% solution non-wetting agent
 - 2.67. 100% Lithium Tetraborate ($\text{Li}_2\text{B}_4\text{O}_7$) flux
3. PROCEDURE:
 - 3.1. Prepare porcelain crucibles by igniting at 950°C to a constant weight. Cool and store in a desiccator to avoid absorption of moisture.
 - 3.2. Weigh 1.0 gram of cement sample accurately to 0.0001 grams into a prepared porcelain crucible. Ignite sample to a constant weight in a muffle furnace at 950°C and cool in a desiccator.
 - 3.3. Calculate the loss on ignition (LOI) using the following formula:
$$\text{LOI} = (A/B) \times 100$$

where: A = weight of sample after ignition;
B = weight of original sample.
 - ~~3.4. Place at least 1.0 gram of cement sample from LOI in the spex 8000 mixer/mill and grind until a fine powder is obtained. This takes approximately 30 seconds.~~

- 3.45. Weigh accurately to 0.0001 grams 6.0 grams of flux directly into platinum crucible. Then weigh accurately to 0.0001 grams 0.6 grams of LOI free cement sample ~~from the Spex 8000 mixer/mill~~ directly into the platinum crucible. Add 3 drops of LiBr solution. Place the platinum crucible and dish in the Perl'X 3 machine and select to run program 9* for all types of cement. This takes approximately 15 minutes. Program run is dependent on sample type.
- 3.65. ~~Open the cement application and e~~Enter ~~names of sample~~the sample identification -and LOI information in ~~the~~ measure sample screen on the measure and ~~—————~~analyze program.
- 3.7-6. Place sample in a 27mm steel cup. Then place in x-ray instrument and prepare to run cement application on the measure sample screen. Click on measure at the bottom of the screen. This may take a few minutes. The application chosen is dependent upon the sample type.
4. QUANTIFICATION: Program quantifies data by using a least squares program. Similar samples with known chemical make-ups are used as standards in the quantification technique. As many standards as possible are used for best quantification. The results are reported as oxides in weight percents.
5. REPORT:
 - 5.1. % LOI
 - 5.2. % Insolubles
 - 5.3. % SO₃
 - 5.4. % Al₂O₃
 - 5.5. % Fe₂O₃
 - 5.6. % C₃A (2.650 x %Al₂O₃) – (1.692 x %Fe₂O₃)
 - 5.7. % MgO
 - 5.8. % Na₂O
 - 5.9. % K₂O
 - 5.10. % CaO
 - 5.11. % SiO₂
 - 5.12. % Total Alkali (Na₂O + 0.658 x K₂O)

* NOTE: Program 9 includes the following: One oxidation for 2 minutes, temperature 1100°C, power of generator 77, agitation angle 25, and agitation speed 10. One fusion for 6 minutes, temperature 1100°C, power of generator 77, agitation angle 50 and agitation speed 15. Then there is a pause before casting for 10 seconds at a temperature of 1100°C. Casting lasts 2 minutes, temperature 1100°C, casting angle 123, casting speed 10 and time for solidification is 30 seconds. Lastly, there is natural cooling for 1 minute and forced air cooling for 3 minutes at a flow rate of 40. The setting of the dish height dial is 40/12 and this depends on the size of the platinum dish used.

APPROVED

DIRECTOR
DIVISION OF MATERIALS

DATE

03/07/08

~~APPROVED~~

~~Director~~

~~DIVISION OF MATERIALS~~

~~DATE~~

~~1/6/03~~

Kentucky Method 64-254-083

Revised ~~1/6/03~~ 03/07/08

Supersedes 64-254-030

Dated ~~1/27/00~~ 01/06/03

Km254083.doc