

## TEST METHOD FOR ANALYSIS OF EPOXY CONTENT

### 1. SCOPE

- 1.1 This method is designed to quantitatively determine the epoxy content of the epoxy component of two component structural steel coatings.

### 2. MATERIALS AND REAGENTS

- 2.1 Balance accurate to 0.001g
- 2.2 Burette, 10 ml capacity with 0.05 ml graduations
- 2.3 Graduated cylinder, 25 ml
- 2.4 Erlenmeyer flasks, 125 ml
- 2.5 Volumetric flasks, 50 ml
- 2.6 Magnetic stirrer, adjustable speed
- 2.7 Magnetic stirring bars, PTFE coated
- 2.8 High speed centrifuge
- 2.9 Centrifuge tubes, 50 ml, disposable polypropylene
- 2.10 Xylene, reagent grade
- 2.11 Chloroform-Chlorobenzene mixture (1:1), reagent grade
- 2.12 Crystal Violet indicator solution
- 2.13 Glacial Acetic Acid, reagent grade
- 2.14 Hydrogen Bromide (HBr), 30% by weight in acetic acid (approx. 5.78 N)
- 2.15 Potassium acid phthalate ( $\text{KHC}_8\text{H}_4\text{O}_4$ ), primary standard grade
- 2.16 Potassium acid phthalate ( $\text{KHC}_8\text{H}_4\text{O}_4$ ) 0.1N solution

### 3. PREPARATION OF STANDARD SOLUTIONS

- 3.1 Crystal Violet Indicator Solution:  
Dissolve 0.1 g crystal violet in 100 ml glacial acetic acid.
- 3.2 Potassium acid phthalate ( $\text{KHC}_8\text{H}_4\text{O}_4$ ) 0.1N solution:  
Dry the potassium acid phthalate in an oven @  $110^\circ\text{C}$  to a constant weight; cool in dessicator. Dissolve 1 g (weighed to 0.001 g) potassium acid phthalate in 25 ml of glacial acetic acid by gently heating. Record weight. Cool to room temperature.

### 4. STANDARDIZATION OF HYDROGEN BROMIDE SOLUTION

- 4.1 The 30% HBr solution should be standardized each day.
- 4.2 To the potassium phthalate solution (see 3.2) add 10 drops of crystal violet indicator solution.
- 4.3 Titrate the potassium acid phthalate solution with the 30 % HBr solution to a blue-green end point.

### 5. SAMPLE PREPARATION: (Perform in duplicate)

- 5.1 To a tared centrifuge tube add 5 g sample, and record weight to 0.001 g.
- 5.2 Thin the sample with 5 ml xylene.
- 5.3 Tap the bottom of the centrifuge tube on the side of the counter to disperse the sample.
- 5.4 Centrifuge the sample to produce a relative centrifugal force (RCF) of 4200 for 30 minutes. See 8.1 for calculation of RCF.
- 5.5 Transfer the supernatant to the 125 ml Erlenmeyer flask in which the titration will be performed.
- 5.6 Re-disperse the pigment in the bottom of the centrifuge tube with 5 ml of xylene.
- 5.7 Centrifuge the sample again for 30 minutes as in 5.4.
- 5.8 Add this supernatant to the original supernatant for analysis.

### 6. PROCEDURE: (Perform in duplicate)

- 6.1 To the supernatant in the 125 ml Erlenmeyer flask add 25 ml of a 1:1

mixture of chloroform and chlorobenzene.

- 6.2 Place a magnetic stirring bar into the flask and mix on the magnetic stirrer to dissolve; stir at a moderate speed to avoid splashing.
- 6.3 Add 10 drops of crystal violet indicator solution.
- 6.4 Titrate with the 30 % HBr in acetic acid solution to a blue-green endpoint with the stirrer rotating at a moderate speed to avoid splashing.
- 6.5 Slow down the titration near the endpoint to allow ample time for the reaction to take place (see 7.1).
- 6.6 Record the amount of titrant used.
- 6.7 Make a blank determination on the reagent in an identical manner.
- 6.8 Record the amount of titrant used for the blank determination.

## 7. NOTES

- 7.1 After each addition allow time for the reaction to take place. The endpoint is a blue-green color. When approaching the endpoint additions should be dropwise. If the solution turns yellow the sample has been over-titrated.

## 8. CALCULATIONS

- 8.1 Calculation for Relative Centrifugal Force (RCF):

$$RCF = 0.00001118 \times r \times N^2$$

Where: RCF = Relative Centrifugal Force  
r = Radius of Rotation  
N = Revolutions Per Minute (RPM)

- 8.2 Calculation for Normality of the 30% HBr in acetic acid:

$$N = (W \times 1000) / (204.2 \times V)$$

Where: W = potassium acid phthalate used, g  
V = HBr solution used, ml

- 8.3 Calculation for weight per epoxy equivalent (grams of that component which contains one gram equivalent of epoxy groups):

$$WPE = 1000W / N(V-B)$$

Where: V = HBr solution used for titration of the sample, ml  
B = HBr solution used to titrate the blank, ml  
N = normality of the HBr solution  
W = weight of the sample used, g

9. REPORT

- 9.1 Report the weights of the duplicate coating samples to the nearest 0.001g.
- 9.2 Report the Normality of the HBr solution.
- 9.3 Report the WPE in duplicate in grams.
- 9.4 Results should be considered suspect if duplicate analyses differ by more than 5% relative percent difference when performed by a single analyst.

APPROVED

\_\_\_\_\_  
DIRECTOR  
DIVISION OF MATERIALS

DATE

\_\_\_\_\_  
03/31/08

Kentucky Method 64-264-08  
Dated 03/31/08

KM26408.doc